

SurgiSuite

Multi-functional surgical platform



**Far Infrared
Animal Warming Pad
with Controller
(optional)**

*Pad recessed in center
of platform for a level
surgical field*

Chemical Resistant Base

*Stainless steel surgical surface
and polycarbonate overlay for
added protection*

**Bright LED Lighting
(optional)**

*Connects to multi-port
USB power strip*

SurgiSuite shown with optional add-ons

Order your SurgiSuite:

sales@apacscientific.com

APAC
SCIENTIFIC



Kent Scientific
CORPORATION

All SurgiSuite includes:

- Easy to clean, chemical-resistant surface
- Magnetic anesthesia mask/intubation tubing stabilizer, (1-small and 1-large)
- Magnetic limb positioners (4)
- Elastomer wire positioner (1-meter)
- Detachable procedure bowl
- Disposable platform covers (5-pack)

1-Year 100% Parts & Labor Warranty

Specifications:

	SURGI-M Mouse	SURGI-R Rat
Base dimensions:	12.90 in x 16.80 in (32.77 cm x 42.67 cm)	15.75 in x 17.80 in (40.01 cm x 45.21 cm)
Stainless steel platform:	11.90 in x 7.90 in (30.23 cm x 20.07 cm)	12.75 in x 10.75 in (32.39 cm x 27.31 cm)
Warming surface:	6.6 in x 5 in (16.51cm x 12.7cm)	8.75 in x 6.75 in (22.23 cm x 17.15 cm)

Configurations:



Surgical board
SURGI-M / SURGI-R



Surgical board with warming
SURGI-M02 / SURGI-R02



Surgical board with LED Lighting and Magnifier Kit
SURGI-M01 / SURGI-R01



Surgical board with warming, LED Lighting and Magnifier Kit
SURGI-M03 / SURGI-R03



Surgical board with warming, LED Lighting and Magnifier Kit and Intubation Kit
SURGI-M04 / SURGI-R04

Both the RightTemp® Jr. temperature monitor & homeothermic control and the PhysioSuite® with RightTemp and other optional physiological monitoring modules can easily connect to the SurgiSuite to provide regulated temperature control.

Order your SurgiSuite today:
sales@apacscientific.com

APAC
SCIENTIFIC



Kent Scientific
CORPORATION

Far Infrared Warming Products

Choose the best solution for your research



Kent Scientific's Far Infrared warming products provide a variety of options to safely warm your animals. From simple pad power regulation to automatic control based on temperature feedback, we have the solution to every research warming requirement.

Features & Benefits:

- **Deep Penetration**
Warms animal quickly and safely
- **Small Footprint**
Easily fits into your surgical setup
- **Homeothermic options**
Regulate animal temperature



*Woven PVC fabric warming pad
available in 2 sizes*

Order your Far Infrared products today:

sales@apacscientific.com

APAC
SCIENTIFIC



Kent Scientific
CORPORATION

Far Infrared Warming Product Selection Guide



Warming Pad Control

- Small footprint; easily fits into your surgical set-up
- Large display; quickly set and identify target temperature
- Operates using AC power

Starting at \$375



RightTemp® Jr Homeothermic Control

- Automatic control measures and regulates your animal's temperature
- Includes temperature sensor
- Far infrared warming pad provided

Starting at \$725



PhysioSuite® RightTemp® Module with Advanced Features

- Temperature monitor and homeothermic control
- Modular physiological monitor system
- Optional modules for pulse oximetry and heart rate and end-tidal CO₂

Starting at \$975

	Warming Pad Control	RightTemp® Jr Homeothermic Control	PhysioSuite®
Homeothermic control	○	●	●
Pad temperature displayed	●	●	●
Body temperature displayed	○	●	●
Temperature sensor(s) included	●	●	●
15.2cm x 20.3cm woven PVC fabric pad	●	●	○
20.3cm x 25.4cm woven PVC fabric pad	●	●	●
Internal Memory	○	○	●
Data export capability	○	○	●
Real-time display	●	●	●
User custom menu	○	○	●
Adjustable temperature settings	○	●	●
AC powered	●	●	●
Electrophysiology compatible	○	○	●
Optional CO ₂ module	○	○	●
Optional pulse oximetry module	○	○	●

Order your Far Infrared products today:

sales@apacscientific.com
apacscientific.com

APAC
SCIENTIFIC



Kent Scientific
CORPORATION